



Laerdal[®]
helping save lives



Bystander CPR Saves Lives

Little Anne[®] QCPR

Providing feedback based on just watching students perform can be challenging for instructors. The next generation of our Little Anne[®] manikin includes QCPR measurement and feedback technology that takes the guesswork out of instruction—all from the palm of your hand. Little Anne QCPR works with an Instructor app that allows you to monitor and coach up to six students at once, a Learner app that lets students monitor their own performance, and a Classroom app for training large groups.

Little Anne QCPR helps instructors improve CPR training quality and efficiency in their courses, while adding a fun and engaging element for students that promotes learning and motivation!



laerdal.com



Little Anne® QCPR

Little Anne® QCPR is still the same durable, lightweight manikin as the original Little Anne, but now includes QCPR technology that provides:

- Real-time feedback on compressions and ventilations
- A summary with an overall score and improvement suggestion
- Details on compression release, depth and rate, ventilation volume, number of compressions/ventilations/cycles

Feedback Options

When used with your choice of feedback option(s), Little Anne QCPR meets the AHA's requirement for feedback devices.



QCPR Instructor app

For instructors to monitor multiple students



- Free download on App Store and Google Play
- Connects to 1-6 manikin(s)
- Shows live and summative feedback of CPR performance
- Provides suggestion for improvement, post-session
- Results are saved in the app
- QCPR Race game reinforces quality CPR skills in a fun and engaging way



QCPR Learner app

For students to monitor their own performance



- Free download on App Store and Google Play
- Connects to 1 manikin
- Shows live and summative feedback of CPR performance
- Provides suggestion for improvement, post-session
- Results are saved in the app



QCPR Classroom app

For instructors to monitor large groups



- Free download on App Store - iPad only
- Connects to 1-42 manikin(s)
- Shows live and summative feedback of CPR performance
- Provides group performance summary
- QCPR Race game reinforces quality CPR skills in a fun and engaging way



SkillGuide™

- Provides a tethered, "plug and play" option
- Can be used standalone or as a supplement to either of the apps

Little Annes manufactured after February 1999 can be quickly and easily upgraded to QCPR with the Little Anne QCPR Upgrade Kit.



Ordering Information

| | |
|-----------|---|
| 123-01050 | Little Anne QCPR |
| 123-03050 | Little Anne QCPR Dark Skin |
| 124-01050 | Little Anne QCPR 4-pack |
| 124-03050 | Little Anne QCPR 4-pack Dark Skin |
| 123-60750 | Little Anne QCPR Upgrade Kit – compatible with Little Annes manufactured after Feb 1999 |

Accessories and Consumables

| | |
|-----------|---|
| 123-30050 | SkillGuide with extension cable |
| 020300 | Little Anne Airways (pkg. 24) |
| 020301 | Little Anne Airways (pkg. 96) |
| 15120103 | Manikin Face Shields (6 rolls of 36 ea) |
| 152400 | Manikin Wipes (pkg. 50) |
| 152401 | Manikin Wipes (pkg. 1200) |
| 310210 | Manikin Faces (pkg. 6) |
| 310215 | Manikin Faces, Dark Skin (pkg. 6) |
| 123-40050 | QCPR Race Circles, 6-Pack |

Visit Laerdal.com/LittleAnneQCPR to learn more about Little Anne QCPR and how you can upgrade your existing Little Anne manikins.

