



For a healthy  
start in life



## 8-CHANNEL HEAD ARRAY COIL 1.5 / 3.0T

- 8-channel Array, cylindrical, receive only coil
- High-resolution MR Imaging of the brains of newborns and premature babies with a head circumference of up to approx. 40cm
- Approved for use inside MR Diagnostics  
Incubator System **nomag**<sup>®</sup> IC **ADVANCED**

Compatible with leading MRI devices like Siemens, Philips, GE  
Further details upon request

Change of technical specification without notice



Kindly provided:  
Oslo University Hospital, Rikshospitalet, Norway

