

**microlife<sup>®</sup>**

A partner for people. For life.

Leading  
Technologies

Smart solutions for  
better health.



Explore the product range of Microlife

**Product catalogue**



<b>Blood Pressure</b>	<b>8</b>
Fully-automatic upper arm blood pressure monitors	12
Fully-automatic wrist blood pressure monitors	19
Semi-automatic upper arm blood pressure monitors	22
Cuffs and supplies	23
Overview blood pressure monitors	24
Manual aneroid blood pressure kits	26
Stethoscopes	28
<b>Blood Glucose Meter</b>	<b>30</b>
Blood glucose monitoring system	34
<b>Fever</b>	<b>42</b>
Bluetooth® thermometers	44
Infrared thermometers	45
Digital thermometers	52
Overview infrared and digital thermometers	58
<b>Respiratory Care</b>	<b>60</b>
Compressor nebulisers	62
Mesh nebulisers	72
Asthma monitor	73
Overview nebulisers	74
Fingertip pulse oximeters	76
<b>Baby Care</b>	<b>80</b>
Breast pumps	82
Overview breast pumps	85
Accessoires	86
Electric nasal aspirator	88
<b>Weight</b>	<b>90</b>
Diagnostic scales	92
Weight scales	95
<b>Microlife Connected Health+</b>	<b>96</b>
<b>Professional Products</b>	<b>98</b>
Office blood pressure monitors	102
Home blood pressure monitors	104
Ambulatory blood pressure monitors	106
Cuffs and supplies	108
Overview professional blood pressure monitors	110

“Our vision is to empower people through our products and services to be able to manage and promote their own health.”



## A partner for people. For life.

Microlife is one of the world leaders in the development and manufacturing of medical diagnostic equipment for home monitoring and institutional use. Blood pressure monitors, blood glucose monitoring system, digital thermometers and respiratory care devices are Microlife's core business.

## Empowering people

Microlife is the world's leading manufacturer of thermometers. The company is also one of the world market leaders in the development and production of medical diagnostic equipment – for use at home and in health facilities. Blood pressure monitors, digital thermometers, respiratory care devices and flexible heat pads are some of our core products.

## Our guiding principle

As the pioneer of the mercury-free digital thermometer, Microlife products are free of environmentally hazardous substances. We strive for high quality, reliability and accuracy. We develop technology for the benefit of your health and set the quality standard just as high as for professional applications.

## Our products

Our medical devices are approved for hospital use and guarantee the highest level of safety. They are robust and user-friendly. Whether it's fever, blood pressure, respiratory care, weight management or flexible heating, our technology offers precise and reliable solutions.

## Our focus: Home diagnostics

We supply everything you need to analyse your condition at home. People want to take responsibility for themselves. We support you in this. This is our future.



### National Institute for Health and Care Excellence (NICE) Recommendation

Microlife WatchBP Home A is the only blood pressure monitor officially recommended BY NICE, for use during routine blood pressure measurement for all GP's and primary care professionals in the United Kingdom (NHS).



reddot

### Red Dot Design Award

A part of our blood pressure monitors and digital thermometers have been recognised with one of the world's most reputable design award, the "red dot" product design.

“In the future, we want to continue working together to do what we do best: Providing strong solutions for people. For life.”



# microlife<sup>®</sup> World Map



**Regional headquarters**  
Taiwan  
Switzerland  
USA

**Branch office**  
China  
Lithuania  
Russia  
Netherlands  
France  
United Kingdom  
USA

**Manufacturing site**  
China  
Vietnam

**microlife<sup>®</sup>**

A partner for people. For life.



# Blood Pressure

---

Knowing your blood pressure and knowing if you have Atrial Fibrillation are the first steps in stroke prevention.



## Blood pressure technologies



### MAM technology

MAM (Microlife Average Mode) is the analysis of 3 successive blood pressure measurements in less than 3 minutes so that superior reliability is achieved by automatic weight data. MAM technology provides highest accuracy for blood pressure monitoring both at home and in professional settings. Multiple readings give important information to doctors and ensures that treatment plans meet healthcare needs.



### AFIBsens technology

Microlife blood pressure monitors with implemented AFIB technology, allow patients to be screened for Atrial Fibrillation during blood pressure measurements at home, using a sophisticated algorithm. If AF is detected during blood pressure measurements, the AFIB symbol is displayed.



### IHB technology

The IHB (Irregular heartbeat) technology analyses pulse frequency and detects irregular heartbeat during blood pressure measurement. If an irregular heartbeat occurs during the measurement, the IHB symbol is displayed. IHB does not replace a cardiac examination, but it serves to detect pulse irregularities at an early stage.



### Gentle+

The Gentle+ technology improves the inflation and deflation rate during blood pressure measurements. It is an adaptive pressure control technology, which allows optimal speed and pressure control thereby ensuring a comfortable measurement.



### Comfort+

Comfort+ technology or also named IMT (Inflation Mode Technology), is a rather new measurement technique and allows an extra gentle, quick and accurate blood pressure measurement already during inflation. Thereby the patient suffers less long from the unpleasant pressure on the upper arm and usually, the measurement duration is shorter.



### Cuff Fit Check

Microlife Cuff Fit Check technology evaluates the signal during cuff inflation of a measurement, to determine if the cuff fit is sufficiently tight or not, and then displays corresponding symbols to indicate the cuff fit status. Please note that Cuff Fit status will be displayed only after the cuff is inflated to pressure above the user's diastolic pressure, because it takes time to check the pulse signal.



### Microlife Connected Health+

is a useful app for easily monitoring your blood pressure and other important health values like temperature or weight.



### BP Analyzer+

This software helps you to get more comprehensive information regarding your blood pressure readings.



### Clinical validation

passed according to the ESH protocol and ISO81060-2:2013.



### Clinical validation

passed according to the BIHS protocol with A/A grade. Additionally, these monitors are validated for use in children from the age of 12 and special patient groups such as pregnancy and pre-eclampsia, renal disease, obesity and diabetes.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Microlife Corp. is under license. Other trademarks and trade names are those of their respective owners.





▶ Product video

## BP A7 TOUCH BT

Bluetooth® blood pressure monitor with stroke risk detection and touch screen



### Stroke risk detection

High accuracy in Atrial Fibrillation (AF) detection for stroke prevention



### Microlife Average Mode (MAM)

3 measurements in 1 click and average result for superior reliability



### App and PC software included

The Connected Health+ app and the BPA+ software enable you to easily analyse your blood pressure over a long-term



reddot

### Features

- AFIBsens technology
- Bluetooth
- PC software BPA+ included
- 2 users with guest mode
- Memory-average function
- Traffic light indicator
- Lock function
- Cuff check indicator
- MAM technology
- Microlife Connected Health+ app included
- Gentle+
- 2 x 99 data memory
- Adjustable time and date
- Rechargeable batter-friendly
- Low battery indication
- Arm movement indicator

[www.microlife.com/a7-touch-bt](http://www.microlife.com/a7-touch-bt)



▶ Product video

# BP B6 Connect



Bluetooth® blood pressure monitor with stroke risk detection



### Stroke risk detection

High accuracy in Atrial Fibrillation (AF) detection for stroke prevention



### Microlife Average Mode (MAM)

3 measurements in 1 click and average result for superior reliability



### App and PC software included

The Connected Health+ app and the BPA+ software enable you to easily analyse your blood pressure over a long-term

### Features

- Fully-automatic for upper arm
- Atrial Fibrillation detection with AFIBsens
- IHB for detection of irregular heartbeat
- Optimal speed and pressure control with Gentle+
- LED traffic light indicator for blood pressure classification
- Extra large display
- Cuff fit check
- Arrow average indicator "MyCheck"
- 2 users
- 2 x 99 data memory
- Average of stored measurements of the last 28 days
- Clinical blood pressure average "MyBP"
- Movement indicator
- Adjustable date and time
- Low battery indicator
- Rechargeable battery-friendly

[www.microlife.com/b6-connect](http://www.microlife.com/b6-connect)



## BP A6 BT



Bluetooth® blood pressure monitor with stroke risk detection



**Stroke risk detection**  
High accuracy in Atrial Fibrillation (AF) detection for stroke prevention



**Microlife Average Mode (MAM)**  
3 measurements in 1 click and average result for superior reliability



**App and PC software included**  
The Connected Health+ app and the BPA+ software enable you to easily analyse your blood pressure over a long-term

### Features

- AFIBsens technology
- MAM technology
- Bluetooth
- Microlife Connected Health+ app
- USB cable and BPA+ software included
- Gentle+
- 99 data memory
- Memory-average function
- Traffic light indicator
- Switch for lock/unlock
- Adjustable date and time
- Low battery indication
- Cuff check indicator
- Arm movement indicator

[www.microlife.com/a6-bt](http://www.microlife.com/a6-bt)



▶ Product video

# BP B3 AFIB



Clinically Validated  
BHS protocol: A/A

Blood pressure monitor  
with stroke risk detection



**Stroke risk detection**

High accuracy in Atrial Fibrillation (AF) detection for stroke prevention



**Microlife Average Mode (MAM)**

3 measurements in 1 click and average result for superior reliability



**Conical soft cuff M-L (22–42 cm)**

Comfortable cuff with washable cuff cover for longlasting, hygienic use

**Features**

- Fully-automatic for upper arm
- Extra large display
- AFIBsens technology
- IHB for detection of irregular heartbeat
- Optimal speed and pressure control with Gentle+
- Cuff fit check
- Traffic light indicator for blood pressure classification
- 99 data memory
- Average of stored measurements
- Adjustable date and time
- Low battery indicator
- Rechargeable battery-friendly

[www.microlife.com/b3-afib](http://www.microlife.com/b3-afib)



▶ Product video

# BP B3 Comfort PC

Advanced blood pressure monitor with Comfort+ for an extra-smooth measurement during cuff inflation



**3 MAM** **Microlife Average Mode (MAM)**  
3 measurements in 1 click and average result for superior reliability

**Comfort+**  
For a quick, accurate and extra comfortable measurement

**PC connection and software**  
PC connection with BP Analyzer+ software and USB cable included

### Features

- Fully-automatic for upper arm
- Comfort+ technology
- Software BPA+ included for PC analysis
- MAM technology
- 2 users with 2 x 99 data memory each
- Average of stored measurements
- Cuff fit check
- Traffic light for easy BP classification
- Adjustable date and time
- Low battery indicator
- Rechargeable battery-friendly

[www.microlife.com/b3-comfort-pc](http://www.microlife.com/b3-comfort-pc)



▶ Product video

# BP B2 BASIC

Blood pressure monitor  
with irregular heartbeat detection



Clinically Validated  
BHS protocol: A/A



**Detection of irregular heartbeat**  
Irregular heartbeat (IHB) detection for early alerts to possible heart irregularities



**Gentle+**  
Optimal speed and pressure control for comfortable measurements



**Washable soft cuff M-L (22–42 cm)**  
Comfortable cuff for all common arm sizes

## Features

- Fully-automatic for upper arm
- IHB technology
- Gentle+ technology
- Cuff fit check
- 30 data memory
- Average of stored measurements
- Easy indicator for blood pressure classification
- Adjustable date and time
- Low battery indicator
- Rechargeable battery-friendly

[www.microlife.com/b2-basic](http://www.microlife.com/b2-basic)



# BP B1 CLASSIC

Blood pressure monitor with optimal speed and pressure control



Clinically Validated  
BHS protocol: A/A



**Detection of irregular heartbeat**  
Irregular heartbeat (IHB) detection for early alerts to possible heart irregularities.



**Gentle+**  
Optimal speed and pressure control for comfortable measurements



**30 data memory and average display**  
Automatic storage of the last 30 blood pressure values and display of mean value of all stored readings

**Features**

- Fully-automatic for upper arm
- Easy indicator for blood pressure classification
- Soft cuff M-L (22-42 cm)
- 30 data memory
- Cuff fit check
- Movement indicator

[www.microlife.com/b1-classic](http://www.microlife.com/b1-classic)



# BP W10 Basic

Wrist blood pressure monitor with irregular heartbeat detection



Clinically Validated  
BHS protocol: A/A



### Detection of irregular heartbeat

Irregular heartbeat (IHB) detection for early alerts to possible heart irregularities.



### Cuff Fit Check

Irregular heartbeat (IHB) detection for early alerts to possible heart irregularities.



### 30 data memory

Automatic storage of the last 30 readings

### Features

- Cuff fit check
- Pre-formed cuff (13,5 – 21,5 cm)
- 1 user
- 30 data memory
- Date / time
- Follow-up model for BP W1 Basic

[www.microlife.com/w10-basic](http://www.microlife.com/w10-basic)



# BP W70 BT

Wrist blood pressure monitor with Comfort+ inflation measurement



Bluetooth®



**Detection of irregular heartbeat**  
Irregular heartbeat (IHB) detection for early alerts to possible heart irregularities.



**Cuff Fit Check**  
Irregular heartbeat (IHB) detection for early alerts to possible heart irregularities.



**99 data memory**  
Automatic storage of the last 99 readings

**Features**

- Inflation measurement with Comfort+
- Silent measurement
- Wrist position check
- Bluetooth® connection
- 2 users
- 99 memories

[www.microlife.com/w70-bt](http://www.microlife.com/w70-bt)



# BP w3 COMFORT

Wrist blood pressure monitor with Comfort+ inflation measurement



Clinically Validated  
ISO protocol



### Detection of irregular heartbeat

Irregular heartbeat (IHB) detection for early alerts to possible heart irregularities.



### Comfort+

Extra gentle measurement during the inflation of the cuff



### Extra-large wrist cuff

For wrist circumferences between 13.5cm up to 23cm

### Features

- Inflation measurement with Comfort+
- Irregular heartbeat (IHB) detection
- 30 data memory and average display
- Soft wrist cuff with D-ring (13.5-23cm)
- BP classification with single-colour eco traffic light
- Adjustable date and time
- Cuff fit check
- Rechargeable battery-friendly

[www.microlife.com/w3-comfort](http://www.microlife.com/w3-comfort)



## BP N2 EASY

Semi-automatic blood pressure monitor with IHB technology



### Detection of irregular heartbeat

Irregular heartbeat (IHB) detection for early alerts to possible heart irregularities.



### One button operation

Easy to use



### Conical soft cuff M (22–32 cm)

Comfortable conical cuff

### Features

- IHB technology
- 1 data memory
- Hypertension warning
- One-button operation
- Rechargeable battery-friendly

[www.microlife.com/n2-easy](http://www.microlife.com/n2-easy)

# Cuffs and supplies

Cuffs, as important as the blood pressure monitor



## Rigid cuff

Conical, preshaped upper arm cuff for high level of comfort

### Features

- Wide-range cuff size
- Conical, fan-shaped cuff for optimal fit
- Compatible with all Microlife blood pressure monitors
- Rigid material for easy self-attachment
- Clinically validated
- Available size:  
M-L (22–42 cm/8.7–16.5")

[www.microlife.com/rigid-cuff](http://www.microlife.com/rigid-cuff)

## Soft cuff

Conical, soft upper arm cuff with washable cuff cover

### Features

- Washable cuff cover for long-term and hygienic use
- Conical, fan-shaped cuff for optimal fit
- Compatible with all Microlife blood pressure monitors
- Metal bar for easy self-attachment
- Clinically validated in all sizes
- Different sizes available:  
S (17–22 cm/6.7–8.7")  
M (22–32 cm/8.7–12.6")  
M-L (22–42 cm/8.7–16.5")  
L (32–42 cm/12.6–16.5")  
L-XL (32–52 cm/12.6–20.5")

[www.microlife.com/soft-cuff](http://www.microlife.com/soft-cuff)

## Overview

Find your blood pressure monitor



**BP A7 Touch BT**



**BP B6 Connect**



**BP A6 BT**



**BP B3 AFIB**

	Data memory	2 × 99	2 × 99	99	99	
	Atrial fibrillation detection	•	•	•	•	
	Microlife Average Mode	•	•	•	•	
	Detection of Irregular Heartbeat		•		•	
	Gentle+	•	•	•	•	
	Comfort+					
	Bluetooth® and app	•	•	•		
	PC connection and software	•	•	•		
	Traffic light indicator	•	•	•	•	
	One button operation					
	2 users	•	•			
	Average display	•	•	•	•	
	Cuff fit check		•		•	
	Adapter included	•	•	•		
	Guarantee	5 years	5 years	5 years	5 years	
	Full Cuff Range (S, M, M-L, L-XL)	•	•	•	•	



BP B3 Comfort PC    BP B2 Basic    BP B1 Classic    BP W10 Basic    BP W70 BT    BP W3 Comfort    BP N2 Easy

	2x99	30	30	30	99	30	1
	•						
	•	•	•	•	•	•	•
		•	•				
	•				•	•	
					•		
	•	•	•	•	•	•	
							•
	•				•		
	•	•	•	•	•	•	
	•	•	•	•	•	•	
	5 years						
	•	•	•				

no power socket



## BP AG1-40

---

Ergonomic and professional aneroid blood pressure kit



### Precise deflation control

Integrated dial and pumping ball provides good support during inflation and deflation



### Professional cuff (25–40 cm)

Soft wide-range cuff for professional fitting



### Ergonomic operation

Large scale dial for easy and professional reading

### Features

- Manual aneroid blood pressure kit
- High-quality sphygmomanometer
- Stethoscope
- Professional cuff (without metal bar)
- Bulb and valve
- Soft bag

[www.microlife.com/ag1-40](http://www.microlife.com/ag1-40)



# BP AG1-20

---

Cost-effective, basic  
aneroid blood pressure kit



**Complete solution**

Stethoscope included in the kit



**Easy to use**

Easy to read, inflate and deflate



**Soft cuff M (22–32 cm)**

Soft cuff for normal fitting

**Features**

- Manual aneroid blood pressure kit
- High-quality sphygmomanometer
- Stethoscope
- Bulb and valve
- Soft bag

[www.microlife.com/ag1-20](http://www.microlife.com/ag1-20)



# ST 77

---

Lightweight Sprague  
Rappaport stethoscope with  
turnable double-chestpiece



High sound precision



Reliable quality



Accessories included

#### Features

- Stethoscope
- 3 spare ear tips of different sizes
- 3 bell type chest pieces
- 2 spare transparent plastic diaphragms
- Available in 5 colours:  
black, dark blue, gray, red, lilac

[www.microlife.com/st-77](http://www.microlife.com/st-77)



# ST 71

---

Single head light-weight stethoscope with flat chest piece

- ✓ High sound precision
- ✓ Reliable quality
- ✓ Lightweight chest piece

#### Features

- Stethoscope
- Lightweight aluminium chest piece
- Adjustable chrome-plated brass binaural
- Durable one-piece moulded PVC "Y" tubing
- Ear tips with soft cone-shaped bell

[www.microlife.com/st-71](http://www.microlife.com/st-71)

# microlife<sup>®</sup>

A partner for people. For life.



A person wearing a green sweater and a brown knit hat is sitting at a wooden table in a rustic kitchen. The background features a brick wall, a window with greenery outside, and a white electric kettle on a wooden countertop. Three mugs are hanging on a rack on the brick wall.

## Blood Glucose

---

High precision, less pain  
and innovation technology  
for monitoring blood  
glucose easily.

## 2 blood glucose technologies you should know



### AC/PC setting

**AC (Ante Cibum)** and **PC (Post Cibum)** meal markers provide more details on the test results of your blood glucose before and after meals. Testing your blood glucose both before and after meals, allows you to see how that meal affects your blood glucose levels and helps you to understand which meals may be best for your blood glucose control.

During the day, levels tend to be at their lowest just before meals. A normal blood sugar level **is less than 100 mg/dL** after not eating (fasting) for at least 8 hours. **And less than 140 mg/dL** 2 hours after eating.

The most powerful influence on blood glucose levels comes from food. Whether you have **type 1 or type 2 diabetes**, the peak blood glucose levels are often likely to occur around two hours after a meal. The advantage of testing in so called "pairs" of before and after meals, is that you can see how much your sugar levels have gone up (or even down) between meals.




### HCT – Hematocrit compensation

**Hematocrit (HCT)** is the Volume percentage of red blood cells in the blood, measured as part of a blood test. The measurement depends on the number and size of red blood cells. The hematocrit is slightly more accurate, as the packed cell volume (PCV) includes small amounts of blood plasma trapped between the red cells. HCT levels are estimated by measuring oxidation of mediator diffusion to the upper gold electrode during the first portion of the assay. The HCT-corrected glucose levels are determined by an on-meter algorithm.

**The technology used in the Microlife Gluco-A+ provides accurate blood glucose measurements that are insensitive to HCT levels across the range 20-60%.**





## Gluco-A<sup>+</sup>

Blood glucose monitoring system with gold electrode



**Hematocrit compensation**  
Increasing accuracy for wide HCT range patients



**AC/PC Settings**  
Before/After meal measurement setting



**Risk indicator**  
Ideal for use at home and on the move

### Features

- For capillary and venous blood testing
- Wide hematocrit measurable range
- Gold electrode
- 4 Alarm settings
- No coding required
- Automatic turn on/off
- Measurement in 7 seconds
- Instant 7, 14, 30, 90 days average
- 1000 results storage

[www.microlife.com/gluco-a](http://www.microlife.com/gluco-a)



# Gluco-A<sup>+</sup> Test Strips

For Microlife Gluco-A<sup>+</sup>



GDH-FAD reagent



7 seconds testing time



1 µL blood sample

## Features

- Microlife Gluco-A<sup>+</sup> Glucose Test Strips are used with the Microlife Gluco-A<sup>+</sup> Glucose Meters.
- They are intended for self-testing by people with diabetes at home and health care professionals in a clinical setting to monitor glucose concentrations in capillary whole blood drawn from the fingertips.
- This test strips provide high accuracy and consistency in testing patients blood glucose, through advanced tri-electrodes glucose sensing technology, with gold electrode test strips.
- Available with 25, 50, 2x25 and 2x50 pieces

[www.microlife.com/glucoa-accessories](http://www.microlife.com/glucoa-accessories)



## GlucoTeq

Blood glucose monitoring system



5  
sec

### Fast measurement

Measurement in 5 seconds



### AC/PC Settings

Before/After meal measurement setting



### 0.5µL blood sample

Requires just a tiny 0.5µL blood sample

### Features

- Strip insert illuminate guidance
- Lighting slot for low ambient luminance
- GDH reagent
- 0.5µL blood sample
- Measurement in 5 seconds
- AC/PC mark
- Automatic turn on/off
- Instant 7, 14, 28, 60, 90 days average
- 500 results storage
- No coding required

[www.microlife.com/glucopeq](http://www.microlife.com/glucopeq)



# GlucoTeq Test Strips

For Microlife GlucoTeq



GDH-FAD reagent

5  
sec

5 seconds testing time



0.5µL

0,5 µL blood sample

## Features

- Microlife GlucoTeq Glucose Test Strips are used with the Microlife GlucoTeq Glucose Meters
- They are intended for self-testing by people with diabetes at home and health care professionals in a clinical setting to monitor glucose concentrations in capillary whole blood drawn from the fingertips
- They are for testing outside the body (in vitro diagnostic use only)
- Available with 25, 50, 2x25 and 2x50 pieces

[www.microlife.com/glucoTEQ-accessories](http://www.microlife.com/glucoTEQ-accessories)



## DiaRite

Blood glucose monitoring system



5  
sec

### Fast measurement

Measurement in 5 seconds

0.5µL

### 0.5µL blood sample

Requires just a tiny 0.5µL blood sample

strip

### Strip ejector

Easy ejection of used test strips

### Features

- Strip ejection
- AC/PC (before/after)
- 0.5µL blood sample
- Measurement in 5 seconds
- 450 results storage
- No coding required
- PC link

[www.microlife.com/diarite](http://www.microlife.com/diarite)



# DiaRite Test Strips

For Microlife DiaRite



GDH-FAD reagent



5 seconds testing time



0.5µL

0,5 µL blood sample

## Features

- Microlife GlucoTeq Glucose Test Strips are used with the Microlife GlucoTeq Glucose Meters.
- They are intended for self-testing by people with diabetes at home and health care professionals in a clinical setting to monitor glucose concentrations in capillary whole blood drawn from the fingertips.
- They are for testing outside the body (in vitro diagnostic use only).
- Available with 25, 50, 2x25 and 2x50 pieces

[www.microlife.com/diarite-accessories](http://www.microlife.com/diarite-accessories)



# Control Solution Level 1

---

For all models of Microlife Blood Glucose Monitoring Systems



For use in self-testing outside the body



For in vitro diagnostic

#### Features

- Microlife Control Solution is for use with all models of Microlife Blood Glucose Monitoring System as a quality control check to verify if the meter is working properly or if your test results are accurate.

[www.microlife.com/gluco-accessories](http://www.microlife.com/gluco-accessories)



# Control Solution Level 3

For all models of Microlife Blood Glucose Monitoring Systems



For use in self-testing outside the body



For in vitro diagnostic

#### Features

- Microlife Control Solution is for use with all models of Microlife Blood Glucose Monitoring System as a quality control check to verify if the meter is working properly or if your test results are accurate.

[www.microlife.com/gluco-accessories](http://www.microlife.com/gluco-accessories)

# *microlife*<sup>®</sup>

A partner for people. For life.



# Fever

---

High precision,  
quality and innovative  
technologies for  
fever management.





## NC 150 BT

Bluetooth® forehead thermometer for easy temperature monitoring



3  
sec

### 3-seconds measurement

Fast and easy measurement without touching the person

### Bluetooth® technology



Instant temperature monitoring on your smartphone

36.2°C  
37.6°C

### Silent Glow™ technology

Green: reading up to 37.4°C  
Red: reading higher than 37.4°C

### Features

- Infrared forehead thermometer
- Enables easy monitoring of the temperature by using Bluetooth® Smart
- Connects to the “Microlife Connected Health+” app (available free for iOS and Android)
- Blue tracking light – easy to use at night
- Multifunctional (body, object and ambient temperature)
- 30 data memory (stored with date and time)
- Sound control – silent mode option for measurements at night
- Clinically tested and recommended by doctors

[www.microlife.com/nc-150-bt](http://www.microlife.com/nc-150-bt)



▶ Product video

# NC 200

Non contact thermometer with auto-measurement and distance control



reddot award 2018 winner



### 3-seconds measurement

Touch-free measurement in only 3 seconds



### Auto measurement with distance control

Takes a measurement automatically, when a distance within 5 cm is detected



### Guidance system for self-measurement

The self-indicator light on the back shows the user that the device is at the right distance

### Features

- Infrared forehead thermometer
- Precise, gentle and hygienic
- 30 data memory (stored with date and time)
- Multifunctional (body, object and ambient temperature)
- Blue tracking light – easy to use at night
- Illuminated fever alarm: Silent Glow™ technology
- Sound control – silent mode option for measurements at night
- Clinically tested in Europe
- Red Dot Award: Product Design 2018

[www.microlife.com/nc-200](http://www.microlife.com/nc-200)



## NC 300

Non contact thermometer  
with an easy to use ergonomic design



3  
sec

### 3-seconds measurement

Fast and easy measurement  
without touching the person



### Blue tracking light

A blue LED light indicates the correct  
measurement area and the distance

36.2°C  
37.6°C

### Silent Glow™ technology

Green: reading up to 37.4 °C  
Red: reading higher than 37.4 °C

### Features

- Infrared forehead thermometer
- Ergonomically shaped design – very easy handling
- Precise, gentle and hygienic
- 30 data memory (stored with date and time)
- Multifunctional (body, object and ambient temperature)
- Sound control – silent mode option for measurements at night
- Clinically tested in Europe

[www.microlife.com/nc-300](http://www.microlife.com/nc-300)



# NC 400

Non contact forehead thermometer with child-friendly design



3  
sec

### 3-seconds measurement

Fast and easy measurement without touching the person



### Blue tracking light

A blue LED light indicates the correct measurement area and the distance



### Silent Glow™ technology

Green: reading up to 37.4 °C  
Red: reading higher than 37.4 °C

### Features

- Infrared forehead thermometer
- Unique dolphin design - specifically developed to make a measurement more pleasant for children
- Precise, gentle and hygienic
- Multifunctional (body, object and ambient temperature)
- 30 data memory (stored with date and time)
- Sound control – silent mode option for measurements at night
- Clinically tested in Europe

[www.microlife.com/nc-400](http://www.microlife.com/nc-400)



# IR 310

Ear thermometer with probe cover and ACCUsens guidance system



1  
sec

**1-second measurement**  
Reliable and accurate measurement of the ear temperature

ACCU  
sens ✓

**ACCUsens guidance system**  
Confirms the right position in the ear with a “good” on the LCD display and a beep



**Touch-free probe cover ejector**  
The hygienic and easy way to remove and discard used probe covers

**Features**

- Infrared ear thermometer with probe cover
- Completely safe and hygienic for use by the whole family
- Probe cover detector ensures the probe cover is well-installed for an accurate measurement
- Illuminated probe and display for easy measurements
- Silent Glow™ technology and fever alarm
- 30 data memory (stored with date/time)
- Sound control – silent mode option for measurements at night
- 20 pieces hygienic probe covers included (BPA and latex free)

[www.microlife.com/ir-310](http://www.microlife.com/ir-310)



▶ Product video

# IR 210

Ear thermometer with probe cover and ACCUsens guidance system



1  
sec

### 1-second measurement

Reliable and accurate measurement of the ear temperature

ACCU  
sens ✓

### ACCUsens guidance system

Confirms the right position in the ear with a “good” on the LCD display and a beep



### Hygienic probe covers

BPA and latex free probe covers ensure a hygienic and safe measurement

### Features

- Infrared ear thermometer with probe cover
- Completely safe and hygienic for use by the whole family
- Illuminated probe and display for easy measurements
- Silent Glow™ technology and fever alarm
- 30 data memory (stored with date/time)
- Multifunctional (body, object and ambient temperature)
- Sound control – silent mode option for measurements at night
- 20 pieces hygienic probe covers included (BPA and latex free)

[www.microlife.com/ir-210](http://www.microlife.com/ir-210)



▶ Product video

## IR 200

Ear thermometer with guidance system for accurate positioning



1  
sec

### 1-second measurement

Reliable and accurate measurement of the ear temperature

ACCU  
sens ✓

### ACCUsens guidance system

Confirms the right position in the ear with a “good” on the LCD display and a beep



### Illuminated probe and “clean me” technology

A blue probe LED light enables you to find the correct ear position in the dark

### Features

- Infrared ear thermometer
- Silent Glow™ technology and fever alarm
- 30 data memory (stored with date/time)
- Multifunctional (body, object and ambient temperature)
- Sound control – silent mode option for measurements at night
- Probe cover free
- Clinically tested in Europe

[www.microlife.com/ir-200](http://www.microlife.com/ir-200)





MT 410



MT 3010



## MT 410/3010

World's first antimicrobial thermometer with copper alloy



### Antimicrobial effectiveness

Copper effectively kills a wide range of microbes



### 1 data memory

See the last measurement result



### Beeper

Signal tone once the measurement is completed

### Features

- Digital thermometer with copper alloy
- The antimicrobial copper is applied on the thermometer's tip and battery cover
- 10 seconds measurement (MT 410)
- 60 seconds measurement (MT 3010)
- Flexible tip/soft touch (only MT 410)
- Large illuminated display (only MT 410)
- Automatic switch-off
- Clinically tested and validated
- Worldwide patented technology

Worldwide patent: WO/2013/064847

[www.microlife.com/antimicrobial-thermometer](http://www.microlife.com/antimicrobial-thermometer)



# MT 850

Digital thermometer with  
3 in 1 technology



8  
sec

**8-seconds measurement**  
Accurate, fast and comfortable

3 in 1

**3 in 1 technology**  
Body mode, water bath mode,  
ambient mode



**Gold Temp™ technology**  
Ideal for people who are  
sensitive to nickel

## Features

- Digital thermometer
- Extra large display
- Silent Glow™ technology
- Beeper and fever alarm
- Flexible tip/soft touch
- 100 % waterproof
- Last measurement memory
- Mercury free
- Automatic switch-off
- Low battery indication

[www.microlife.com/mt-850](http://www.microlife.com/mt-850)



## MT 808

Digital thermometer for a rapid and comfortable measurement



**8-seconds measurement**  
Accurate, fast and comfortable



**Flexible measuring tip**  
For a gentle and comfortable measurement



**Fever alarm**  
Beep sound if the reading is higher than 37.4 °C

### Features

- Digital thermometer
- Large display for easy usage
- Beeper
- 100 % waterproof
- Last measurement memory
- Mercury free
- Automatic switch-off
- Low battery indication

[www.microlife.com/mt-808](http://www.microlife.com/mt-808)



# MT 800

Digital thermometer for a fast, easy and comfortable measurement



10  
sec

**10-seconds measurement**  
Accurate, fast and comfortable



**Flexible measuring tip**  
For a gentle and comfortable measurement



**Fever alarm**  
Beep sound if the reading is higher than 37.4 °C

### Features

- Digital thermometer
- Large display for easy usage
- Beeper
- 100 % waterproof
- Last measurement memory
- Mercury free
- Automatic switch-off
- Low battery indication

[www.microlife.com/mt-800](http://www.microlife.com/mt-800)



MT 700



MT 710



MT 720



## MT 700/710/720

Digital thermometer with child-friendly design and flexible tip



30  
sec

**30-seconds measurement**  
Accurate, fast and comfortable



**Ideal for children**  
Appealing design especially made for children



**Flexible measuring tip**  
For a gentle and comfortable measurement

### Features

- Digital thermometer
- Last measurement memory
- Mercury free
- Beeper
- Easy battery change
- Low battery indication
- Automatic switch-off
- Clinically tested

[www.microlife.com/mt-700](http://www.microlife.com/mt-700)

[www.microlife.com/mt-710](http://www.microlife.com/mt-710)

[www.microlife.com/mt-720](http://www.microlife.com/mt-720)



# MT 600

Digital thermometer for a reliable and easy measurement



60  
sec

**60-seconds measurement**  
Accurate and reliable measurement



**100% waterproof**  
Safe and hygienic



**1 data memory**  
See the last measurement result

## Features

- Digital thermometer
- Easy battery change
- Automatic switch-off
- Low battery indication
- Mercury free
- Clinically tested

[www.microlife.com/mt-600](http://www.microlife.com/mt-600)

## Overview

Find your infrared fever thermometer



**NC 150 BT    NC 200    NC 300    NC 400    IR 310    IR 210    IR 200**

Measurement time	3 sec.	3 sec.	3 sec.	3 sec.	1 sec.	1 sec.	1 sec.
Measurement range	0.1–99.9°C	0.1–99.9°C	0.1–99.9°C	0.1–99.9°C	32.0–43.0°C	0.1–99.9°C	0.1–99.9°C
Non contact measurement	•	•	•	•			
Ear measurement					•	•	•
<b>3in1</b> Multifunctional	•	•	•	•		•	•
Bluetooth® and app	•						
Auto-measurement		•					
ACCUsens technology					•	•	•
Blue tracking light	•	•	•	•			
Silent Glow™ technology	•	•	•	•	•	•	•
Silent mode option	•	•	•	•	•	•	•
Illuminated probe					•	•	•
Probe cover					•	•	
Probe cover detector					•		
Probe cover ejector					•		
Illuminated display	•	•	•	•	•	•	•
30 data memory	•	•	•	•	•	•	•
°C/°F switchable	•	•	•	•	•	•	•
Clinically tested	•	•	•	•	•	•	•
5 year guarantee	5 years	5 years	5 years				

## Find your digital fever thermometer



MT 850

MT 808  
MT 800

MT 410

MT 700  
MT 710  
MT 720

MT 3010

MT 600

 Measurement time	8 sec.	8 / 10 sec. (MT 808 / MT 800)	10 sec.	30 sec.	60 sec.	60 sec.
 Measurement range	34.0–43.9°C	32.0–43.9°C	32.0–42.9°C	32.0–42.9°C	32.0–43.9°C	32.0–43.9°C
 Flexible tip	•	•	•	•		
 Fever alarm	•	•	•			
<b>3in1</b> Multifunctional	•					
 Beeper	•	•	•	•	•	•
 Illuminated display	•		•			
 Waterproof	•	•				•
 Gold Temp™ measuring tip	•					
 Antimicrobial properties			•		•	
 Silent Glow™ technology	•					
 Automatic switch-off	•	•	•	•	•	•
<b>0% Mercury</b> Mercury free	•	•	•	•	•	•
 Clinically tested	•	•	•	•	•	•
 Lifetime guarantee	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime

*microlife*<sup>®</sup>

A partner for people. For life.



# Respiratory Care

---

Inhale effectively  
and breathe easy  
with Microlife.





## NEB 310

Professional compressor  
nebuliser 3 in 1



### Efficient treatment

High performance system for the upper and medium/lower respiratory tract

3 in 1

### 3 nebulisation rates

Choose from 3 different nebulisation rates



### Air flow technology

Adjustable flow rate for individual treatment

### Features

- Powerful piston pump for continuous use
- Nasal washer for the treatment of upper airways
- Breath-enhanced nebuliser for lower airways
- Easy one button operation
- Suitable for all common inhalation solutions
- Hygienic and safe – accessories are washable and disinfectable
- Fully equipped with an extra soft adult and child face mask, mouthpiece with expiratory valve, nosepiece and spare filters

[www.microlife.com/neb-310](http://www.microlife.com/neb-310)



▶ Product video

# NEB PRO

Professional compressor nebuliser 2 in 1



### Efficient treatment

High performance system for the upper and medium/lower respiratory tract

2 in 1

### Efficient 2 in 1 aerosol therapy

Nasal washer for the treatment of upper airways, breath-enhanced nebuliser for lower airways



### High speed nebuliser

For short inhalation times

### Features

- Powerful piston pump for continuous use
- Easy one button operation
- Suitable for all common inhalation solutions
- Large accessory compartment to keep everything safe and clean
- Cable storage
- Hygienic and safe – accessories are washable and disinfectable
- Fully equipped with an extra soft adult and child face mask, mouthpiece with expiratory valve, nosepiece and spare filters

[www.microlife.com/neb-pro](http://www.microlife.com/neb-pro)

# microlife®



## NEB 150 MINI

Portable compressor nebuliser



### Efficient treatment

For the upper and lower respiratory tract



### One for all

Suitable for all common inhalation solutions



### Small and lightweight

Portable design –  
Ideal for travel

### Features

- Micro-compressor technology
- Suitable for all common inhalation solutions
- Accessories are washable and disinfectable
- Nebuliser holder – keep the nebuliser in its holder on the back of the device, so it will be handy for the next treatment
- Fully equipped with adult and child face mask, mouth- and nosepiece, spare filters and carrying bag

[www.microlife.com/neb-150-mini](http://www.microlife.com/neb-150-mini)



# NEB 210

Compressor nebuliser with air flow technology



### Efficient treatment

For the upper and lower respiratory tract



### Fast treatment

High speed nebuliser for short inhalation times



### Air flow technology

Adjustable flow rate for individual treatment

### Features

- Innovative nebulisation system for efficient treatment
- Suitable for all common inhalation solutions
- Accessories are washable and disinfectable
- Easy storage and transportation with carrying bag
- Nebuliser holder – keep the nebuliser in its holder on the back of the device, so it will be handy for the next treatment
- Fully equipped with adult and child face mask, mouth- and nosepiece, spare filters and carrying bag

[www.microlife.com/neb-210](http://www.microlife.com/neb-210)

# microlife®



▶ Product video

## NEB 200

Compressor nebuliser



### Efficient treatment

For the upper and lower respiratory tract



### Fast treatment

High speed nebuliser for short inhalation times



### Easy one button operation

Suitable for the whole family

### Features

- Innovative nebulisation system for efficient treatment
- Suitable for all common inhalation solutions
- Accessories are washable and disinfectable
- Easy storage and transportation with carrying bag
- Nebuliser holder – keep the nebuliser in its holder on the back of the device, so it will be handy for the next treatment
- Fully equipped with adult and child face mask, mouth- and nosepiece, spare filters and carrying bag

[www.microlife.com/neb-200](http://www.microlife.com/neb-200)



MicroLife®



## NEB NANO

Rechargeable and portable  
compressor nebuliser



### Efficient treatment

For the upper and lower  
respiratory tract



### Rechargeable by micro USB

Battery operated and easily  
rechargeable using USB cable



### Ultra light and silent

For an efficient and comfortable  
on the go treatment

### Features

- Patented micro-compressor technology
- Compact and lightweight design – ideal for travel
- Easy one button operation with battery indicator
- Suitable for all common inhalation solutions
- Easy storage and transportation with carrying bag
- Nebuliser holder for easy use
- Equipped with adult and child face mask, mouthpiece, nosepiece, spare filters, power adapter (medical type) with micro-USB cable and carrying bag

[www.microlife.com/neb-nano](http://www.microlife.com/neb-nano)



# NEB NANO BASIC

Portable compressor nebuliser



### Efficient nebulisation system

For the upper and lower respiratory tract



### Portable design - Ideal for travel

Small, lightweight and ultra-silent - only 180g



### One for all

Suitable for all common inhalation solutions

### Features

- Patented micro-compressor technology
- Easy one button operation
- Accessories are washable and disinfectable
- Easy storage and transportation with carrying bag
- Nebuliser holder for easy use
- Equipped with adult and child face mask, mouthpiece, nosepiece, spare filters, power adapter (medical type) with micro-USB cable and carrying bag

[www.microlife.com/neb-nano-basic](http://www.microlife.com/neb-nano-basic)



## NEB 400

Children's compressor nebuliser



### Efficient treatment

For the upper and lower respiratory tract



### High speed nebuliser

For short inhalation times



### Child-friendly design

For a comfortable and silent operation

### Features

- Fast and comfortable treatment for the whole family
- Suitable for all common inhalation solutions
- One button design – easy use
- Accessories are washable and disinfectable
- Easy storage and transportation with carrying bag
- Nebuliser holder – keep the nebuliser in its holder on the back of the device, so it will be handy for the next treatment
- Fully equipped with adult and child face mask, mouth- and nosepiece, spare filters and carrying bag

[www.microlife.com/neb-400](http://www.microlife.com/neb-400)



# NEB 410

Children's compressor nebuliser



**Efficient treatment**  
For the upper and lower respiratory tract



**High speed nebuliser**  
For short inhalation times



**Child-friendly design**  
For a comfortable and silent operation

## Features

- Fast and comfortable treatment for the whole family
- Suitable for all common inhalation solutions
- One button design – easy use
- Accessories are washable and disinfectable
- Easy storage and transportation with carrying bag
- Nebuliser holder – keep the nebuliser in its holder on the back of the device, so it will be handy for the next treatment
- Fully equipped with adult and child face mask, mouth- and nosepiece, spare filters and carrying bag

[www.microlife.com/neb-410](http://www.microlife.com/neb-410)



## NEB 800

Mesh nebuliser



**Innovative nebulisation system**  
Stainless steel mesh



**Ultra light and silent**  
Ideal for the entire family



**Compact design**  
Battery-operated for comfortable treatment even on the go

### Features

- Portable, pocket size design
- Low medicine residue
- Silent with low power consumption
- One button – user friendly interface
- Supplied with adult and child face masks, mouthpiece, carrying bag and AC adapter (optional)

[www.microlife.com/neb-800](http://www.microlife.com/neb-800)



# PF 100

Asthma monitor



**High accuracy**

Peak expiratory flow (PEF) & FEV<sub>1</sub> measurement



**240 data memory**

Stored with date and time



**Traffic light system**

For an immediate assessment of peak expiratory flow

**Features**

- High accuracy
- Rotating vane wheel method
- Washable flow sensor
- Easy to read LCD display
- Suitable for children and adults
- Asthma analyser software – useful for the management of the patient's respiratory diary
- Supplied with software, connection cable and professional carrying bag

[www.microlife.com/pf-100](http://www.microlife.com/pf-100)

## Overview

Find your nebuliser



**NEB 310**

**NEB PRO**

**NEB NANO**

**NEB NANO BASIC**

**NEB 150 MINI**

<b>Nebulisation rate</b>	mode 1: 0.2/0.3 ml/min. mode 2: 0.3/0.4 ml/min. mode 3: 0.3/0.5 ml/min.	0.37 ml/min.	0.25 ml/min.	0.25 ml/min.	0.3 ml/min.
<b>Particle size</b>	mode 1: ≈ 3.0 μm (MMAD) mode 2: ≈ 4.0 μm (MMAD) mode 3: ≈ 5.0 μm (MMAD)	3.5 μm (MMAD)	2.9 μm (MMAD)	2.9 μm (MMAD)	2.9 μm (MMAD)
<b>Nebuliser capacity</b>	min. 2 ml; max. 10 ml	min. 2 ml; max. 8 ml	min. 2 ml; max. 6 ml	min. 2 ml; max. 6 ml	min. 2 ml; max. 6 ml
<b>Sound noise level</b>	51.5 dBA	52 dBA	45 dBA	45 dBA	52 dBA
<b>Operating cycle</b>	Continuous	Continuous	30 min. on/30 min. off	30 min. on/30 min. off	30 min. on/30 min. off
<b>Weight</b>	1,200 g	1,765 g	200 g	180 g	850 g
<b>Size</b>	160 x 161 x 90 mm	289 x 180 x 127 mm	117 x 69 x 46 mm	117 x 69 x 46 mm	111 x 77 x 165 mm
<b>Guarantee</b>	5 years	5 years	2 years	2 years	5 years
Efficient treatment	•	•	•	•	•
Fast treatment	•	•			
<b>2 in 1</b> 2 in 1		•			
<b>3 in 1</b> 3 nebulisation rates	•				
Adjustable air flow	•	•			
Portable			•	•	•
Rechargeable			•		
Lightweight and silent			•	•	•
Accessories washable	•	•	•	•	•



**NEB 400**



**NEB 410**



**NEB 210**



**NEB 200**



**NEB 800**

	0.4 ml/min.	0.4 ml/min.	0.4 ml/min.	0.4 ml/min.	0.2 ml/min.
	2.83 µm (MMAD)	2.83 µm (MMAD)	≤ 2.08 µm (MMAD)	2.83 µm (MMAD)	≈ 5 µm (MMAD)
	min. 2 ml; max. 8 ml	min. 2 ml; max. 8 ml	min. 2 ml; max. 6 ml	min. 2 ml; max. 8 ml	min. 2 ml; max. 8 ml
	52 dBA	52 dBA	51.5 dBA	52 dBA	<30 dBA
	30 min. on/30 min. off				
	1,424 g	1,300 g	1,200 g	1,306 g	130 g
	200 × 183 × 153 mm	196 × 140 × 134 mm	160 × 161 × 90 mm	103 × 160 × 140 mm	121 × 44 × 59 mm
	5 years	5 years	5 years	5 years	2 years
	●	●	●	●	●
	●	●	●	●	
			●		
					●
					●
	●	●	●	●	●



## OXY 500 BT

Bluetooth® fingertip pulse oximeter



### Medical grade accuracy

Determines oxygen saturation (SpO<sub>2</sub>) and heart rate (pulse)



### Gentle usage

Simple and completely pain-free measurement



### Bluetooth® technology

Instant monitoring on your smartphone

### Features

- Determines oxygen saturation (SpO<sub>2</sub>) and heart rate (pulse)
- Bluetooth® technology for easy monitoring
- Instant monitoring on your smartphone with the Microlife Connected Health+ app
- Dual color OLED display
- Easy to read - Display with 6 available views
- Adjustable display brightness
- Risk reminder
- Low battery indication
- Automatic switch-off
- Supplied with carrying strap and pouch

[www.microlife.com/oxy-500-bt](http://www.microlife.com/oxy-500-bt)



# OXY 210

Fingertip pulse oximeter



**Medical grade accuracy**

Determines oxygen saturation (SpO<sub>2</sub>) and heart rate (pulse)



**Gentle usage**

Simple and completely pain-free measurement



**Risk reminder**

Limit values individually adjustable

**Features**

- Determines oxygen saturation (SpO<sub>2</sub>) and heart rate (pulse)
- Dual color OLED display
- 6 different display modes
- 10 adjustable brightness levels
- Low battery indication
- Automatic switch-off
- Supplied with carrying strap and pouch

[www.microlife.com/oxy-210](http://www.microlife.com/oxy-210)



## OXY 200

Fingertip pulse oximeter



### Medical grade accuracy

Determines oxygen saturation (SpO<sub>2</sub>) and heart rate (pulse)



### Gentle usage

Simple and completely pain-free measurement



### Lightweight and portable

Ideal for use at home and on the move

### Features

- Determines oxygen saturation (SpO<sub>2</sub>) and heart rate (pulse)
- LED display
- Low battery indication
- Automatic switch-off
- Supplied with carrying strap and pouch

[www.microlife.com/oxy-200](http://www.microlife.com/oxy-200)



*microlife*<sup>®</sup>

A partner for people. For life.



# Baby Care

---

Highest comfort with  
Microlife breast pumps.





## BC 100 Soft

Manual breast pump



### Light, easy and ultra-quiet

Allows comfortable expressing of milk anywhere



### BPA free

The bottle is completely free of BPA



### Ideal for newborns

Milk bottle and silicone teat are specially designed for newborns

### Features

- Vacuum technology simulates natural drinking for gentle and comfortable use
- Extra comfort with 2 expressing settings
- Hygienic and safe – accessories are washable and disinfectable
- Fully equipped with 180ml milk bottle (completely BPA free) with cap, bottle stand, two matching screw plugs for bottle and a naturally shaped
- Silicone teat for easy feeding
- Easy storage and transportation with carrying bag

[www.microlife.com/bc100-soft](http://www.microlife.com/bc100-soft)



# BC 200 Comfy

Electric Breast Pump



**Individual usage**

Highest comfort with 10 different stimulation settings



**10 different pumping modes**

Helps mothers to select the most comfortable pump setting for their milk flow



**Easy to use**

Large illuminated display with timer display easily tracks stimulation and expressing time

**Features**

- Designed to mimic baby’s natural sucking rhythm with 2-phase vacuum technology
- Lightweight and compact design for easy expressing at home or on-the-go
- Memory function for personalised usage
- Battery and mains power operated (with battery level indicator)
- Equipped with 180ml milk bottle (completely BPA free) with cap, bottle stand, two matching screw plugs for bottle, naturally shaped silicone teat for easy feeding, medical power adapter and carrying bag

[www.microlife.com/bc200-comfy](http://www.microlife.com/bc200-comfy)



## BC 300 Maxi 2 in 1

Dual electric breast pump



**2in1** Innovative dual mode technology  
Suitable for one and two-sided expressing

**Gentle, convenient and time-saving**  
Double expressing system for more milk in less time

**M** **Highly comfortable**  
Memory function for personalised usage

### Features

- Time-saving with dual mode vacuum technology
- Highest comfort with 10 different stimulation and expressing settings
- Large illuminated display with timer display which easily tracks stimulation and expressing time
- Battery and mains power operated (with battery level indicator)
- Equipped with 2 x 180ml milk bottles with caps, bottle stands, two matching screw plugs for bottle, naturally shaped silicone teat for easy feeding, medical power adapter and carrying bag

[www.microlife.com/bc300-maxi-2in1](http://www.microlife.com/bc300-maxi-2in1)

# Overview

Which breast pump is best for you?



BC 100 Soft



BC 200 Comfy



BC 300 Maxi 2 in 1

Manual breast pump	●		
Single electric breast pump		●	
Dual electric breast pump			●
2-phase vacuum technology	●	●	●
Usage	Occasionally	Often	Frequent
Weight	135 g	330 g	330 g
Size	200 x 160 x 66 mm	130 x 120 x 57 mm	130 x 120 x 57 mm
Guarantee	2 years	2 years	2 years
Stimulation mode	2	10	10
Pumping mode	2	10	10
Memory function		●	●
Display backlight		●	●
Battery and mains power operation		●	●
Light, compact and quiet	●	●	●
BPA free	●	●	●
Washable accessories	●	●	●
Accessories are ideal for newborns	●	●	●



## BC bottle kit

Breast milk bottle



### Compact, light and portable

Allows you to use anywhere



### BPA free

The bottle is completely free of BPA



### Ideal for newborns

Milk bottle and silicone teat are specially designed for newborns

### Features

- 180 ml breast milk bottle with closure
- Ideal to collect, store, freeze and feed breast milk
- Simple to use and clean
- Safe materials for you and your baby
- Dishwasher and microwave safe
- Equipped with 1 x 180 ml breast milk bottle with closure
- Suitable for Microlife BC 100 / 200 / 300

[www.microlife.com/baby-care-accessories](http://www.microlife.com/baby-care-accessories)



# BC replacement & upgrade kit

Suitable for Microlife BC 200 / 300



**Compact, light and portable**

Allows you to use anywhere



**BPA free**

The bottle is completely free of BPA



**Ideal for newborns**

Milk bottle and silicone teat are specially designed for newborns

**Features**

- Vacuum technology simulates natural drinking for gentle and comfortable use
- Safe, easy and comfortable
- Suitable for one and two-sided expressing
- Equipped with 1 x 180ml breast milk bottle with closure cap, 1 x soft silicone cushion, 1 x silicone tube single type, 1 x silicone tube double type
- Suitable for Microlife BC 200 / 300

[www.microlife.com/baby-care-accessories](http://www.microlife.com/baby-care-accessories)

# microlife®



## BC 50

### Electric Nasal Aspirator



0+  
months

#### Ideal for newborns

Suitable for infants from 0 months



#### Waterproof

Easy auto-cleaning function of the pump



#### BPA free

Included soft silicone tips are biocompatible

#### Features

- Gently and quickly removes excess mucus for alleviation of infection of the upper respiratory tract
- 2 comfortable sizes of silicone tips (small and large – suitable for infants and toddlers)
- Battery-operated ergonomic design for comfortable usage
- Compact and lightweight – Ideal for travel
- Quiet and comfortable
- Easy one button operation
- Safe and hygienic – mucus cup and tips are washable and dishwasher-safe

[www.microlife.com/bc-50](http://www.microlife.com/bc-50)



# microlife<sup>®</sup>

A partner for people. For life.



# Weight

---

Our precise weight and diagnostic scales help you to find your healthy weight and keep it.





## WS 200 BT

Bluetooth® diagnostic scale



### Body composition diagnosis

Reliable measurement of your body composition



### Automatic transfer via Bluetooth®

Allows fast analysis of your body composition with the Microlife Connected Health+ app



### ITO coating

ITO (= Indium Tin Oxide) electrode with safety glass weighing platform

### Features

- Analysis of the body weight, fat and water percentage, muscle and bone mass, BMI, with active and basic metabolic rate
- Large, easy to read display
- Automatic transfer via Bluetooth® technology
- Instant weight tracking on your smartphone with "Microlife Connected Health+" app (available free for iOS and Android)
- 180kg capacity
- 100g division
- Kg/St selection

[www.microlife.com/ws-200-bt](http://www.microlife.com/ws-200-bt)



# WS 100

Body fat scale



**Multifunctional**

Measures weight, body fat and body water



**Instant on**

Turn the scale on by tapping it firmly with your foot



**Kg/lb/St selection**

Choose between different weight units

**Features**

- Reliable measurement of your body composition
- Automatically calculates body fat and water percentage in function of height, age and sex
- 10 person memory
- Large, easy to read display
- Athlete mode with 3 levels of activity
- 150kg capacity
- 100g division

[www.microlife.com/ws-100](http://www.microlife.com/ws-100)



## WS 80-N

Diagnostic scale



### Body composition diagnosis

Reliable measurement of your body composition



### Instant on

Step on the scale to begin measurement



Auto off

### Automatic switch-off

Automatic switch-off after a few seconds after the measurement

### Features

- Measures weight, body fat, body water, muscle mass, bone mass, BMI and basal metabolic rate
- Safety glass platform
- Large display
- Slim and compact design
- Easy to use
- 10 person memory
- 2 to 180kg capacity
- 100g division
- Kg/lb/St selection

[www.microlife.com/ws-80-n](http://www.microlife.com/ws-80-n)



# WS 50

Weight scale



#### Instant on

Turn the scale on by tapping it firmly with your foot



Auto off

#### Automatic switch-off

Automatic switch-off after a few seconds after the measurement



#### Kg/lb/St selection

Choose between different weight units

#### Features

- Weight measurement
- Accurate and easy to use
- Slim and compact design
- Safety glass platform
- 150kg capacity
- 100g division

[www.microlife.com/ws-50](http://www.microlife.com/ws-50)



**BP A7 Touch BT**  
[www.microlife.com/a7-touch-bt](http://www.microlife.com/a7-touch-bt)



**BP A6 BT**  
[www.microlife.com/a6-bt](http://www.microlife.com/a6-bt)



**WS 200 BT**  
[www.microlife.com/ws-200-bt](http://www.microlife.com/ws-200-bt)



## Microlife Connected Health+

Helps to check and analyse your health using a few easy steps

**Countless different apps on your smartphone to track your most important health values?**

Who doesn't know that! With the Microlife Connected Health+ you can easily transfer your blood pressure, body temperature and body composition values into one app! This will help to keep you motivated to control them into a positive direction over the long term.



# Microlife Connected Health+

Track your health with Microlife Connected Health+ towards better health.



BP B6 Connect  
[www.microlife.com/b6-connect](http://www.microlife.com/b6-connect)



NC 150 BT  
[www.microlife.com/nc-150-bt](http://www.microlife.com/nc-150-bt)




BP W70 BT  
[www.microlife.com/bp-w70-bt](http://www.microlife.com/bp-w70-bt)



OXY 500 BT  
[www.microlife.com/oxy-500-bt](http://www.microlife.com/oxy-500-bt)

## Data security. Transparent and clear.

Besides your health values, Microlife Connected Health+ also takes the security of your data into account. That's why you decide yourself if you want to use Microlife Connected Health+ locally only, without registration or cloud data storage. If you decide instead to store your data online as well, all relevant information regarding data protection can easily be found directly in the app or on our website – transparent and clear.



# microlife<sup>®</sup>

A partner for people. For life.



# Professional Products

---

Be part of the growing  
movement to prevent  
stroke and heart  
disease.



## WatchBP technologies



### Atrial Fibrillation (AF) detection

AFIB reliably detects Atrial Fibrillation (AF) during blood pressure measurement



### Non-invasive central blood pressure measurement

Non-invasive, reliable method based on brachial pulse volume plethysmography (PVP) waveforms



### Ankle Brachial Index (ABI) assessment

For an easy, fast and accurate assessment of the Ankle Brachial (arm) Index



### Night-time measurement

- For measuring blood pressure during sleep.
- Helps to collect accurate day and night blood pressure measurements for a good overview of the daily blood pressure pattern.



### Double arm measurement

The most reliable tool for determining inter-arm blood pressure differences (IAD)



### Auscultatory mode

For use in e.g. elderly, obese and those with arrhythmia.



### WatchBP average mode

Three consecutive measurements in a row



### Flexible BPM schedule

To help physicians all over the world adhering to measurement recommendations of different hypertension guidelines, Microlife has developed a professional blood pressure monitor with a flexible measurement schedule. With this unique option the professional is free to select the preferred number of measurements, the interval time between measurement, the count-down period and even the way the average blood pressure value is calculated.



### SPRINT algorithm

Performing automated office blood pressure measurements (AOBPM)



### Customisable measurement functions

The device can be upgraded with additional, helpful features like Atrial Fibrillation (AF) detection or central BP measurement function



### PC connection

User-friendly software included



### ESH/AHA/BIHS guidelines embedded

The 7-days-measurement program



### For recording medication intake

The pill-button can be used to record medication intake or to record an "event" that could influence blood pressure level.



### Brachial-ankle pulse wave velocity (baPWV) assessment

User-friendly with high reproducible cuff-based baPWV measurement method to evaluate arterial stiffness. The result is strongly associated with reference carotid-femoral pulse wave velocity (cfPWV), risk factors for arterial stiffness and the presence of arterial damage.





## WatchBP

# Office Vascular

Professional Office blood pressure and cardiovascular screening monitor



Bluetooth®



### Customisable measurement functions

Customise your AOBPM according to your needs with features such as central measurement function or Atrial Fibrillation detection



AFIB

### Atrial Fibrillation (AF) detection\*

AFIBsens detects AF with high accuracy (sensitivity 98%, specificity 92%)



### Non-invasive central blood pressure measurement\*

Non-invasive, reliable method based on brachial pulse volume plethysmography (PVP) waveforms

### Features

- Fast, easy and accurate Ankle Brachial Index (ABI) assessment\*
- Pulse wave velocity (PWV) screening\*
- Dual-cuff design for simultaneous measurement on both arms to assess the Inter-Arm-Difference (IAD).
- Automated office blood pressure measurement (AOBPM) with implemented SPRINT algorithm
- Clinically validated for children aged 3+, pregnancy and pre-eclampsia, diabetes, renal disease and obesity

[www.microlife.com/office-vascular](http://www.microlife.com/office-vascular)

\*Optional features



▶ Product video

## WatchBP Office 2G

Professional automated  
office blood pressure monitor



Bluetooth®



### Customisable measurement functions

Customise your AOBPM according to your needs with features such as central measurement function or Atrial Fibrillation detection



### Atrial Fibrillation (AF) detection\*

AFIBsens detects AF with high accuracy (sensitivity 98%, specificity 92%)



### Non-invasive central blood pressure measurement\*

Non-invasive, reliable method based on brachial pulse volume plethysmography (PVP) waveforms

### Features

- Display of Mean Arterial Pressure (MAP) and Pulse Pressure (PP)
- User-friendly software included
- Easy data transfer via Bluetooth® or USB cable
- Automated office blood pressure measurement (AOBPM) with implemented SPRINT algorithm
- Rechargeable battery pack
- Clinically validated for children aged 3+, pregnancy and pre-eclampsia, diabetes, renal disease and obesity

[www.microlife.com/office-2g](http://www.microlife.com/office-2g)

\*customisable

# microlife®



▶ Product video



## WatchBP

# HOME A

Home blood pressure monitor  
with Atrial Fibrillation (AF) detection



### Atrial Fibrillation (AF) detection

AFIB detects AF with high accuracy  
(sensitivity 98 %, specificity 92 %)



### Usual mode

In usual mode, 3 consecutive measurements  
are automatically performed for high accuracy



### Diagnostic mode

7-day self-measurement program that  
strictly follows ESH/AHA/BIHS guidelines

### Features

- 250 data memory
- PC connection and software included
- Officially recommended by the National Institute for Health and Care Excellence (NICE), UK
- Clinically validated for children aged 12+, pregnancy and pre-eclampsia, diabetes, renal disease and obesity

[www.microlife.com/home-a](http://www.microlife.com/home-a)



▶ Product video

## WatchBP HOME

Home blood pressure monitor



### Usual Mode

For single measurements at any time (250 memory).



### Diagnostic Mode

7-days self-measurement program that strictly follows ESH / AHA / BHS guidelines.



### PC connectivity

Transmit measurement data to PC via USB connection.

### Features

- Home blood pressure monitor with atrial fibrillation (AF) detection function.
- Measuring blood pressure and screening for atrial fibrillation at the same time.
- Easy, reliable and affordable automated blood pressure measurement device.


[www.microlife.com/home](http://www.microlife.com/home)



## WatchBP

# O3 Ambulatory 2G

Professional 24-hour  
blood pressure monitor



Bluetooth®



### Ambulatory blood pressure measurement

User-friendly ambulatory blood pressure monitor for 24-hours blood pressure measurements



### Customisable measurement functions

Customise your ABPM according to your needs with features such as central measurement function or Atrial Fibrillation detection



AFIB

### Atrial Fibrillation (AF) detection\*

AFIB detects AF with high accuracy (sensitivity 98 %, specificity 92 %)

### Features

- PC connectivity with Bluetooth® and USB
- User friendly software for easy data analysis included
- 330 data memory for pill intake, event and blood pressure measurement recordings
- Sophisticated design to fulfil hygiene standards
- Fully-programmable
- Silent mode option
- Clinically validated for children aged 3+, pregnancy and pre-eclampsia, diabetes, renal disease and obesity

[www.microlife.com/ambulatory](http://www.microlife.com/ambulatory)

\*customisable



## Cuffs and supplies

Cuffs, as important as the blood pressure monitor



WatchBP  
**Rigid cuff**

Rigid and conical upper arm cuff, applicable for WatchBP Office devices

**Properties**

- Preformed cuff for highest level of comfort
- All cuff bladders are both latex and PVC free
- Suitable for all common arm sizes
- Available size: M-L (22-42 cm)

[www.microlife.com/watchbp-accessories](http://www.microlife.com/watchbp-accessories)



WatchBP  
**Office cuff**

Soft upper arm cuff, applicable for all WatchBP Office devices

**Properties**

- To fulfil hygienic needs in the professional setting the cuffs are washable in washing machine after removing the bladder
- All cuff bladders are both latex and PVC free
- Available in nylon and cotton
- Available sizes: S (14-22 cm), M (22-32 cm), L (32-42 cm), L-XL (32-52 cm)

[www.microlife.com/watchbp-accessories](http://www.microlife.com/watchbp-accessories)

“Using the wrong size cuff is a major cause for erroneous blood pressure measurement.”



**WatchBP**  
Office ankle cuff

Conical, soft ankle cuff, applicable for WatchBP Office ABI and Central

**Properties**

- Ankle cuff for measuring the Ankle Brachial (arm) Index (ABI)
- All cuff bladders are both latex and PVC free
- 170 cm tube with cuff connector
- The conical shape provides a better fit on special group patients
- Available sizes: M (22–32 cm) and L (32–42 cm)

[www.microlife.com/watchbp-accessories](http://www.microlife.com/watchbp-accessories)



**WatchBP**  
Ambulatory cuff

Soft upper arm ambulatory cuff, applicable for all WatchBP O3 devices

**Properties**

- To fulfil hygienic needs in the professional setting the cuffs are washable in washing machine after removing the bladder
- All cuff bladders are both latex and PVC free
- Available in nylon and cotton
- Available sizes: S (14w–22 cm), M (22–32 cm), L (32–42 cm), L–XL (32–52 cm)

[www.microlife.com/watchbp-accessories](http://www.microlife.com/watchbp-accessories)



**WatchBP**  
Home cuff

Soft upper arm cuff, applicable for all WatchBP Home devices

**Properties**

- With metal bar for easy self-attachment
- All cuff bladders are both latex and PVC free
- Available sizes: S (17–22 cm), M (22–32 cm), M–L (22–42 cm), L (32–42 cm), L–XL (32–52 cm)

[www.microlife.com/watchbp-accessories](http://www.microlife.com/watchbp-accessories)

## Overview

Find your professional blood pressure monitor



WatchBP Office 2G



WatchBP Office Vascular

	Data memory	for 1 patient	for 1 patient
	Atrial Fibrillation detection	optional	optional
	Central BP measurement	optional	optional
	ABI assessment		optional
	Double arm measurement		•
	3 consecutive measurements	•	optional
	SPRINT algorithm implemented	•	•
	24 h ambulatory BP monitor with night-measurement		
	Customisable features	•	•
	PC connectivity and software	•	•
	Pill-button		
	For use in ESRD	•	•
	For use in pregnancy	•	•
	For use in children age 3+	•	•
	For use in children age 12+		
	For use in diabetes	•	•
	baPWV assessment		optional
	Bluetooth® connectivity	•	•



WatchBP Home A



WatchBP Home



WatchBP O3  
Ambulatory 2G

	250	250	
	•		optional
	•		
			•
			•
	•	•	•
			•
	•	•	•
	•	•	•
			•
	•	•	•
	•	•	•
			•

# *microlife*<sup>®</sup>

A partner for people. For life.

4G\_Company\_Brochure\_EN\_4622

Microlife AG  
Esenstrasse 139  
9443 Widnau  
Switzerland



[www.microlife.com](http://www.microlife.com)